

Date d'édition : 12.09.2026

Ref : EWTSOCLA800/4

CLA 800/4 Kit de désinfection à base d'eau, 200 ml

FIELD OF APPLICATION

VERIFORTE® is a rapid and effective water based disinfectant for use on surfaces. VERIFORTE® is pH neutral, alcoholaldehyde free, fragrance and non flammable.

INSTRUCTIONS FOR USE

VERIFORTE® is a ready to use, fluid, clear solution. Surfaces and objects are simultaneously and thoroughly moistened with 50ml/m² using where appropriate a disposable and leaving to dry. Rinsing with water is not necessary. The product evaporates leaving no residue.

HIGHLY EFFECTIVE

VERIFORTE® is tested according to European Standards EN 13697, EN 13704 and EN 14476 and is listed by DGHM/VAH. VERIFORTE® can be used both as a wipe and spray disinfectant. The stated exposure times below, demonstrate the full effectiveness of Veriforte against bacteria, viruses, spores and yeast, even when, mechanical help is not employed.

EXPOSURE TIMES undiluted without mechanical help:

Bacterial Effect surface precleaned 5min*

MRSA??...5 min E.coli???????1 min

E.hirae??..1 min P.aeruginosa??..1 min

Yeasticidal Effect precleaned surface

C. albicans??5 min

*shortest possible exposure time according to VAH - Applied Hygiene Network

Sporocidal Effect precleaned surface

Bac. Subtilis??5 min

Virucidal Effect Murine Norovirus low burden 1 min / high burden 1 min

Adenovirus low burden 1 min / high burden 1 min

Poliovirus low burden 1 min / high burden 30 min

Virucidal against enveloped viruses????..5 min according to the DVV and RKI guidelines (incl. HIV, HBV und HCV)

SPECIFICATIONS Color/Form: transparent/liquid pH value: 6.5-

COMPOSITION >99.5% Water; <0.4% NaCl (Sodium Chloride); <0.05% resp. 0.5g/L NaOCl (Sodium Hypochlorite) / HOCl (Hypochlorous acid)

REGISTRATION N - 54547 for Product type 2 (disinfectants)

N - 54587 for Product type 4 (Food and feed sector)

N - 54588 for Product type 5 (drinking water)