Date d'édition : 01.11.2026

Ref : 688145

Cuvette ronde pour 666 265

Round cuvette for experiments with the UV-C lamp 666266 and the UV-C sensor 666246.
Ozone can be generated in the glass cuvette by applying a high voltage to two opposing needle electrodes. If necessary, the cuvette can be flushed with air via two glass tube extensions that can be closed with plastic caps.

SYSTEMES DIDACTIQUES s.a.r.l.
Savoie Hexapole - Actipole 3 - 242 Rue Maurice Herzog - F 73420 VIVIERS DU LAC
Tel : 04 56 42 80 70 | Fax : 04 56 42 80 71
www.systemes-didactiques.fr